


06/10/2014



Si detectors:  $\sim 1 \text{ mm}$  thick

Bias: -120V

( channel #1  
MPOD slot 6 )

9:

Cable #1: Silicon

Cable #2: Proportional Counter

Cable #3: Master Trigger

Run 1 file: C12-test-beam-001001.root

start: 12:48:42

stop: 13:47:40

Vacuum, no gas. Only getting 10%-20% of the beam

Should only see signals along the beam axis.

Scintillator shows about 250-300 events/s

Run 2 file: C12-test-beam-002001.root

start: 14:19:59

stop: 14:34:32

Optimized beam rate, now get more than scintillator read out

Still no gas in chamber and under vacuum.

Gain for Run 1 & 2 was set to 3.

Run 3 file: C12-test-beam-003001.root

start: 15:40:49

stop: 15:47:53

Chamber filled with methane at 50 torr.

Observed ionization chamber's signal is low, maybe all noise.

Possible causes: gas is not right, voltage is not high enough

Still running at uncalibrated scale due to gain set at 3

instead of 7.

Run 4

file: c12-test-beam-004001.root

start: 15:48:39

stop: 15:54:07

Changed gain to 7. Gas is still at 50 torr.

Run 5

file: c12-test-beam-005001.root

start: 16:35:45

stop: 16:46:41

Pressure changed to 80 torr. gain is set to 7.  
Ionization chamber looks better. No events,  
gas is stopping the beam. Increased intensity of beam  
by factor of 10 so events are seen in silicon. (~~2 events/s~~)  
In the middle of the run, we increased the intensity of  
the beam by another factor of 10. Bias on proportional counter  
was set to 650V.

Run 6

file: c12-test-beam-006001.root

start: 17:09:38

stop: 17:20:21

Pressure change to 70 torr. Proportional counter now on.

Run 7.

file: c12-test-beam-007001.root

start: 17:21:27

stop: 17:31:49

Continuation of run 6.

Run 8

file: c12-test-beam-008001.root

start: 17:53:20

stop: 18:00:55

Pressure at 50 torr

Run 9

file: c12-test-beam-009001.root

start: 18:47:41

stop: 18:49:44

Pressure at 400 torr. Closed shutter near beginning of run.  
Never Ended shortly after closing shutter.

Run 10

file: c12-test-beam-010001.root

start: 19:49:44

stop:

Run 11

file: c12-test-beam-011001.root

start: 19:53:44

stop: 20:09:17

Continuation of run 10.

Run 12

file: c12-test-beam-012001.root

start: ~~19~~: 20:10:52

stop: 20:28:30

Changed gain of proportional counters from 0 to 5.

Stop run, happy ☺

06/11/2014

Determined that either the detectors are not centered or the chamber  
\* is tilted.

Run 13 File: c12-test-beam - 013001.root

Stat: 17:38:17

Stop:

pressure: 350 torr

USE For beam  
E determination

Discovered that the proportional counters  
do not function when the bragg peak  
of the incoming ions occurs within the  
cell. With 360 torr, signals from the  
prop counter were small and extended.  
Reducing pressure to 350 torr the signals  
increased significantly in the proportional  
counter in size.

\* DO NOT STOP BEAM IN PROP \*  
COUNTER

Run 14 File: c12-test-beam - 014001.root

Stat: 18:06:52

pressure: 400 torr

Run 15 File: c12-test-beam - 015001.root

Continuation of the last run with  
intensity ~~roughly~~ 10x higher.

390 stops in proportional counter  
--- changing to 395 torr.

## Current Settings for Module 1 (Si)

Gain: 15, 7, 7, 7  
Thresholds: 255 for 0, 20 for all else  
Shaping: 3, 3, 3, 3,

For Module 0 (prop) we are running in common mode

Gain = 5  
Shaping = 2

→ PRODUCTION FINALLY!!

### Run 16

Start: 20:05:00

File: cl2-test-beam-016001.root

pressure: 397 torr

Moving down to counting room

Grigory

850 - 491 - 6377

Gerveny

850 - 345 - 7587

### Run 17

File: cl2-test-beam-017001.root

pressure: 357 torr

~~DAQ~~ DAQ Crashed!!!

### Run 18

File: cl2-test-beam-018001.root

pressure: 397 torr



\* Run 18 is junk \*  
Some problems recovering  
day.

Run 19

Start: 21:06:16

pressure: 397 torr / 405 torr

File: c12-test-beam-019001.root

SEEMS OK AFTER RESET

Steps to start a run:

or not if  
want to accumulate

① IN Analysis Manager, make sure  
that All Spectra are cleared with  
the "clear arrays" button. on  
the triumph analysis tab under histograms.

② Press "Start AQQ" in the Transport  
manager

a) If writing to disk, make sure  
run number has incremented  
and file prefix is correct

b) supervisor must be entered

③ Make sure that "begin" is pressed  
on analysis manager.

To stop run :

- ① Press stop ~~key~~ on transport manager

To start And program:

- ① Navigate to `~/triumf-backend`
- ② type And
- ③ Open And only after launching TMR !

To start Tmr program:

- ① ssh into `grogaupe@cyfegro1`
- ② navigate to `~/triumf-frontend`
- ③ type Tmr usb

Calculating peak energies:

From run 13 w/ 350 torr  $^{12}\text{C}$  peak appears at  $\sim 2596$  keV

$\Rightarrow 78.497$  MeV beam energy

$\Rightarrow 74.033$  after window

Resonance energies CM Frame:

① 421.4 keV  $J = \frac{1}{2}^+$

② 1558.49  $J = \frac{3}{2}^-$

③ 1603.49  $J = \frac{5}{2}^+$



Penetration distance:

①	<del>48.18</del> cm	42.48 cm
②	<del>42.43</del> cm	37.41 cm
③	<del>42.14</del> cm	37.16 cm

Run 20

file: c12 - test\_beam - 00001.root  
Start: 23:05:52  
press.: 397 torr / 405 torr

back to calculation

Proton	scattering energies	After loss
①	1.556 MeV	<del>1.543 MeV</del>
②	5.754 MeV	0.873 MeV
③	5.92 MeV	5.336 MeV
		5.5015 MeV

High energy peak in 2.3 @ 5.4 MeV

Low energy " " @ 0.9 MeV

CONSISTANT!

If peaking beam, replace  
attenuators ← for AT

Ask to peak if  $I_{cr} < 4000$

Sthan: 574-807-1848

Run 21

file: c12-test-beam-021001.root

start: 00:42:39

press: 397 torr / 405 torr

Run 22

file: c12-test-beam-022001.root

start: 02:16:57

press: 397 torr / 405 torr

Run 23

file: c12-test-beam-023001.root

start: 03:57:25

press: 397 torr / 405 torr

Run 24

file: c12-test-beam-024001.root

start: 05:38:24

press: 397 torr / 405 torr

Run 25

file: c12-test-beam-025001.root

start: 07:18:44

press: 397 torr / 405 torr

stop: 08:51:36

✓ Increased gas pressure to 1020 torr, switch ion chamber out of trigger

### Run 26

Start: 11:21:51

pressure: 1020 torr / 1000 torr

File: c12-test-beam-024001.root

LAST RUN ~~ION~~ PROP COUNTER HAD WRONG SHAPING TIME. STOPPED.

### Run 27

Start: 11:29:48

pressure: 1020 / 1000

File: c12-test-beam-024001.root

GAIN SEEMS TO HIGIT STOPPED + changed from 3 to 5

### Run 28

JUNK

### Run 29

Reset gain back to 5

pressure: 1020/1000

start: ~~11:29:18~~ 12:29:10

File: c-12-test-beam-028001.root

Seems that gain is too low, increasing U band

Run 30

increased bias ~~from~~ by 200V  
Started run but gain was  
to low. [Reran ~~from~~ with  
gain set to 3 instead of 5.]



Run 31

file: c12\_test\_beam - 030001.root  
pressure: 1020/1000  
Start: ~~31~~ 12:52:24  
Stop: 13:16:32

Down Scale factor during test run  
= 1000